

Jasco J-810 CD – Quick ON/OFF Guide

Turning ON the system:

1. Open **Nitrogen Controller** (on Desktop) select the correct COM port (COM'X':N2 valve). Turn the N₂ ON, check that the flow indicator is showing N₂ flow.
2. **Purge the system for at least 10 min** (40 min if it hasn't been used in over 2 days); steps 3 & 4 can be done while waiting.
3. **Turn on the water bath and the temperature controller**; press **Start** on the temperature controller front panel.
4. **Switch on the spectrometer** (once the system is initialised the shutter LED will light up, proceed to step 5 only after this happened)
5. Open **Spectra Manager** (on Desktop), select the type of experiment you want to do.
6. **Turn on** the Xe lamp. On the instrument, the lamp indicator turns green.
7. **The lamp will require at least 10 min to stabilize after ignition**
Note: If it is necessary to restart the lamp after running, it must be allowed to cool for about 15 mins before re-ignition is attempted. **Do not switch off the lamp if there are less than 8 hours between you and the next user.** Contact the next user via e-mail.

Turning OFF the system:

1. **Turn off the Xe lamp.** Make sure the lamp indicator on the instrument has switched off
2. **Quit the software** (there may be two/three windows open. Close all of them before going to step 3)
3. **Switch off the spectrometer**
4. **Switch off the temperature controller and water bath**
5. **Let the system be purged with nitrogen for 10 min.**